

## BX51 Start-up

- 1) Fluorescence only, turn on the mercury burner
- 2) Turn on the camera
- 3) Confirm the camera filter is in the desired position (i.e. grayscale or RGB)
- 4) Turn on the microscope
- 5) Turn on the computer
- 6) Start the camera software
- 7) Confirm the DIC analyzer is out of the lightpath
- 8) For transmitted light imaging, establish Köhler illumination

## BX51 Objectives

All are

U = Universal (Brightfield, Darkfield, DIC, and Polarized Light)

PlanFI = plan-semi-apochromat = corrected for 4 colors chromatically and spherically

Objectives	NA	Resolution $R=0.61*\lambda/NA$ obj
UPlanFL 4x/0.13	0.13	2.581 $\mu\text{m}$
UPlanFL 10x/0.30	0.3	1.118 $\mu\text{m}$
UPlanFL 20x/0.50	0.5	0.671 $\mu\text{m}$
UPlanFL 40x/0.75	0.75	0.447 $\mu\text{m}$
UPlanFL 100x/1.30 oil/iris	1.3	0.258 $\mu\text{m}$

## Shut-down – only shut down if there are no users scheduled within an hour

- 1) Turn off the camera software
- 2) Shut-down the computer
- 3) Turn off the camera
- 4) Set the camera filter to grayscale
- 5) Remove the DIC analyzer from the lightpath
- 6) Turn off the microscope
- 7) Remove oil from the 100x if used.
- 8) Position the 4x or 10x objective over the microscope stage.
- 9) Always turn the mercury burner off last
- 10) Cover the microscope with the dust cover avoiding the mercury lamp housing

## Pixel size by Objective for Images using the Retiga 2000R QImaging Camera at 1600 x 1200

4x = 541 pixels = 1000.00 microns = 1.85 microns per pixel

10x = 1352 pixels = 1000.00 microns = 0.74 microns per pixel

20x = 1702 pixels = 630 microns = .37 microns per pixel

40x = 1892 pixels = 350 microns = .185 microns per pixel

100x = 1351 pixels = 100 microns = .074 microns per pixel