



## BX60 Start-up

- 1) Turn on the mercury burner if you are using it.
- 2) Turn on the camera
- 3) Turn on the microscope
- 4) Establish Köhler illumination
- 5) Turn on the computer and camera software



## BX60 Objectives

All are

U = Universal (Brightfield, Darkfield, DIC, and Polarized Light)

PlanFI = plan-semi-apochromat = corrected for 4 colors chromatically and spherically

Objectives	NA	Resolution $R=0.61 \cdot \lambda / NA \text{ obj}$
UPlanFL 4x/0.13	0.13	2.581 $\mu\text{m}$
UPlanFL 10x/0.30	0.3	1.118 $\mu\text{m}$
UPlanFL 20x/0.50	0.5	0.671 $\mu\text{m}$
UPlanFL 40x/0.75	0.75	0.447 $\mu\text{m}$
UPlanFL 100x/1.30	1.3	0.258 $\mu\text{m}$

## Shut-down – only shut down if there are no users scheduled within an hour

- 1) Turn off the camera software and computer
- 2) Turn off the camera
- 3) Turn off the microscope
- 4) Position the 4x over the microscope
- 5) Always turn the mercury burner off last
- 6) Cover the microscope with the dust cover avoiding the mercury lamp housing

## Pixel size by Objective for Images using Macrofire camera at 2048 x 2048

4x = 552 pixels = 1000.00 microns = 1.81 microns per pixel  
10x = 1387 pixels = 1000.00 microns = .72 microns per pixel  
20x = 1744 pixels = 630 microns = .36 microns per pixel  
40x = 1936 pixels = 350 microns = .18 microns per pixel