

Olympus IX81 Inverted

Epifluorescence & Brightfield
Automated Multi-D image acquisition

Introduction to the
NRI-MCDB Microscopy Facility
IX81 Inverted Microscope

Contents

- Start-up
- Preparing for Imaging
 - Part 1 – General
 - Part 2 – Transmitted – Brightfield
 - Part 3 – Fluorescence
 - Part 4 – Multi-D camera acquisition
- Shut-down




Step 1: Sign-in

- Record the following:
 - Date:
 - Your name:
 - Your Project Code (i.e. Index Code):
 - Your Principal Investigator (PI):
 - Extension (Optional):
 - Time-in (the time you arrived):
 - Time-off (the time you left):
 - Comments (any notes on system condition)




Step 2: Turn on the Microscope / Computer

- A: Sola Light Engine
 - Switch on back
 - Used for fluorescence
- B: Olympus IX2-UCB
 - Automated scope controls
- C: LG-PS2
 - Brightfield light source
- Turn-on the computer
- Log-in using your ADS account name and password.
- For access to the network drive, select Run and then type
- `\\microscopy-nas1.nri.ucsb.edu`
- Create a shortcut for future use.
- Open Micro-Manager 1.4
- Choose configuration file “Olympuslx81Ben.cfg” unless you’ve made your own.




Preparing for Viewing and Imaging

- Part 1 – General Preparation
- Part 2 – Transmitted Light Applications
- Part 3 – Fluorescence Applications
- Part 4 – Multi-D Acquisitions





General Prep – Initial Settings

- The camera, filters, objectives, shutters, light path, and z-position can all be controlled through Micro-Manager.
- Light source changes automatically with FilterCube/Shutter selection.
- There's a footswitch that controls the fluorescence shutter






Part 2: Preparing the microscope for transmitted imaging

- Direct light to the transmitted path.
 - Set FilterCube to “Position-1”
 - Set LightPath to “side port” for camera or “eyepiece”
- Adjust the light intensity.
- If “Auto shutter” is enabled the shutter will open and close depending on the selection of “live” camera
- Disable “Auto shutter” to allow direct control of the shutter
- For single images, use “Snap” or “Album” for more complex operations use “Multi-D Acquisition”




Part 2: Establish Kohler Illumination

- Place a slide on the stage
- Rotate the 10x objective into position.
- Use the microscope focus knob to bring the sample into crisp focus.
- Rotate the field diaphragm to the closed position.
- Use the condenser focus knob to adjust the condenser height so that the field diaphragm appears crisply focused when viewed through the microscope.
- Use the centering knobs located on the left and right sides of the condenser to center the view of the field diaphragm.
- open the field diaphragm just beyond the field of view
- Kohler illumination is objective specific




Part 3: fluorescence imaging

- Choose desired FilterCube (DAPI, CFP, FITC, TRITC)
- Shutter automatically opens and closes with live viewing
- Footswitch will also open and close Spectra shutter
- Adjust SpectraLightLevel to desired brightness to control photobleaching.
- Adjust exposure time to achieve appropriate brightness



Part 4: Multi-D image acquisition

- Time points: settings for a time series
- Z-stacks: if “relative Z” is chosen then positions are relative to current position
- Channels: lets you cycle between different scope configurations.
“FilterCube/Shutter” lets you choose filter settings and light source
- “Acquisition order” is important



Shut-Down Procedure

- Check the online schedule
 - Shut-down if nobody is scheduled within the next hour
 - Leave the system on if somebody is using the system in the next hour but do the following.
 - Log-off the computer
 - Close the fluorescent shutter
 - Clean-up
 - Return to the 10x objective
 - Sign-off in the log.
- Adjust your online reservation end-time if you finished early
- Shut off the computer
- Turn off A, B, and C
- Complete the paper log by filling-in
 - Time you finished
 - Any comments
- Put dust cover over microscope

Specifications

- 4 objectives
 - 10x/0.4 UPlanSApo
 - 20x/0.7 UPlanApo
 - 40x/0.6 LUCPlanFLN
 - long working distance, coverslip thickness adjustment collar
 - 100x/1.3 oil UPlanFL
 - oil immersion, NA adjustment collar
 - 4 fluorescence filters
 - DAPI
 - CFP
 - FITC
 - TRITC
 - Qimaging Scientific CMOS, 12-bit, 1280x1024, 30 fps full resolution
 - Model #01-ROL-BOLT-M-12
- Camera Calibration
- 0.725 $\mu\text{m}/\text{pixel}$
- 0.361 $\mu\text{m}/\text{pixel}$
- 0.178 $\mu\text{m}/\text{pixel}$
- 0.0722 $\mu\text{m}/\text{pixel}$

